

# Isothermal Titration Calorimetry (ITC):

## Theory and Applications

## Instrumentation Basics

**Marco Coletti**

*Sr. Manager – EMEA Scientific Operations*

Waters™



# ITC Theory

# Basics of an Isothermal Titration Calorimetry (ITC)


**Calorimetry** is an application of the First Law of Thermodynamics to heat transfers and allow us to **measure** the **enthalpy** of a reaction including **binding**. It is a **label free** and **non-optical** method meaning no specific molecular groups or labels are required.

- The signal in a calorimeter is the “heat” of a process  $\Rightarrow \Delta H$

$\Delta H$  reflects the amount of heat energy required to achieve a particular state – *direct measurement*


- **Gibbs Free Energy Change** ( $\Delta G$ ): it dictates the direction of binding equilibria;  $\Delta G = -RT \ln K_a$  – *direct measurement*
- We can also determine the value of  $\Delta S$  because it also measures the value of  $K_a$  or  $\Delta G$ .

$\Delta S$  measure how easily this energy is distributed among various molecular energy levels – *indirect measurement*



**Gibbs Free Energy** ( $\Delta G$ ) is a balance between **Enthalpy** ( $\Delta H$ ) and **Entropy** ( $\Delta S$ )


$$\Delta G = \Delta H - T\Delta S$$




**Association Constant,  $K_a$**  vs **Dissociation Constant,  $K_d$** :  $K_d = 1 / K_a$

# Basics of an Isothermal Titration Calorimetry (ITC) Experiment


- Mix two solutions
  - Mixing chemicals will generate heat
- Measure Heat (q)
  - Proper controls ensure certain heat events are measured
- Fit data to assumed model
  - Model based from shape of curve and/or background knowledge
- Analyze data
  - Quantify  $K_d$ ,  $\Delta G$ ,  $\Delta H$ ,  $\Delta S$ ,  $n$ , CMC,  $\Delta C_p$ , enzyme kinetics



Instrument applies or withdraws power to keep  $\Delta T$  constant between reference and sample cell




Raw output is the change in power ( $dP$ ,  $\mu W$ ) with time






- As titrant is introduced to the titrand, instrument will apply or withdraw power
  - Will keep reference and sample cells at constant temperature
- Areas of the peaks are the heat ( $q$ ), and at constant pressure  $q = \Delta H$

# Isothermal Titration Calorimetry

- The titrant (ligand) is loaded into a syringe at 10x greater concentration than the titrand (substrate) in the cell.



- Substrate concentration is ideally in the range of  $1-1000 \cdot K_d$ .
- Heat ( $q$ ) is produced when the complex is formed. As free substrate  decreases, so does the heat .



Thermodynamic parameters:

Stoichiometry,  $n$     Binding constant,  $K_a$ ,    Enthalpy,  $\Delta H$

$$\Delta G = -RT \ln(K_a) \quad \text{and} \quad \Delta G = \Delta H - T\Delta S$$

Free energy,  $\Delta G$     Entropy,  $\Delta S$

# Beyond $K_d$ – Values from ITC

## Quantify:

$K_d$ ,  $\Delta H$ ,  $n$  (stoichiometry)

\*Binding constants in the order of  $10^2$ - $10^9$   $M^{-1}$  can be accurately calculated

$\Delta G$ ,  $\Delta S$

$\Delta C_p$ ,  $\Delta[H^+]$

CMC


## Rationalize:

Structural changes

Lead optimization

- Thermodynamic parameters
  - $\Delta G = -RT \ln(K)$
  - $\Delta G = \Delta H - T\Delta S$
- Heat capacity
  - $\Delta C_p = (\partial\Delta H/\partial T)_p$
- Enzyme kinetics
  - $K_m$ ,  $K_{cat}$
- Other experiments
  - CMC – critical micelle concentration
  - Dimer dissociation

# Consider the contributions to binding



## ■ Enthalpy

- Directly associated with the strength and number of interactions broken or formed
  - Non-covalent interactions
    - Hydrogen bonds, electrostatic interactions, van der Waals forces
- Solvent contributions are important


## ■ Entropy

- Hydrophobic interactions
- Solvation effects (release of water molecules)

## Enthalpy-Entropy Compensation

$$\Delta G = \Delta H - T\Delta S$$

$$\Delta G = -RT\ln(K_a)$$



Waters™



# Instrument

## Nano ITC



## Affinity ITC



## Affinity ITC Auto






## Fixed Cylindrical Cell



## Specifications

- Operating temperature: 2 to 80 °C
- Minimum Detectable Heat: 0.04  $\mu$ J
- Maximum Measurable Heat: 5,000  $\mu$ J
- Low Noise Level: 0.0013  $\mu$ W
- Baseline Stability: 0.02  $\mu$ W/hr
- Cell volume: 190  $\mu$ L (*also* in 1.0 mL)
- Cell geometry: Fixed cylindrical
- Cell material: 24 K Gold (*also* in Hastelloy)
- Temperature Stability: 8  $\mu$ °C at 25 °C
- Temperature Control: Active heating & cooling
- Response Time: 3.3 Seconds



- **FlexSpin Hardware:**
  - *Thin, flexible* support wire
  - **Innovative stir paddle** design
  - **Stirring mechanism separate** from *injection system*
  - Easy removal and replacement

- **FlexSpin Control:**
  - Easy, intuitive on/off and stir speed settings thru ITCRun software
  - Stir speed range – 0-200 rpm
  - Recommended stir speed =125 rpm


- **FlexSpin:**
  - Sharper peaks
  - Faster return to baseline
  - Eliminates dead zones
  - Tempered steel

Waters™



# ITC Applications

- Binding Interactions
- Direct or Competition
- Dimer Dissociation
- Micelle or liposome characterization
- Quality Control, [active]
- Michaelis-Menten Enzyme Kinetics



# Application: Binding Interactions

**Purpose:** Measure the binding affinity between a protein and antibody

## Method:




- 32  $\mu\text{M}$  Protein was titrated into 8  $\mu\text{M}$  IgG (black)
- 32  $\mu\text{M}$  Protein was titrated into buffer (red)

## Observations:

- The ITC technology generates a complete thermodynamic profile of the interaction
- Strong binding affinity ( $K_d = 72 \text{ nM}$ )
- Enthalpy and stoichiometry provide additional information

*ITC is a powerful analytical tool for characterizing the fundamental driving forces of molecular binding reactions*

## Protein-Protein Interactions



The shape of the binding curve determines the accuracy of  $K_a$  and  $\Delta H$

# Application: Binding Interactions

**Purpose:** Measure the binding affinity between a DNA and small molecule

## Method:

- 450  $\mu\text{M}$  Doxorubicin was titrated into 250  $\mu\text{M}$  (bp) DNA Oligo (Blue)
- 450  $\mu\text{M}$  Doxorubicin was titrated into buffer (Red)

## Observations:

- Doxorubicin is one of the most important anti-cancer chemotherapeutic drugs, being widely used for the treatment of solid tumors and acute leukemias
- Average binding affinity ( $K_d = 1.2 \mu\text{M}$ )
- Molar ratio of 0.2 provides insight to how many base pairs the drug binds to

*ITC is a powerful analytical tool for characterizing the fundamental driving forces of molecular binding reactions*

## Small Molecule-DNA Interactions





Image from: <https://doi.org/10.1016/j.bbcan.2013.12.002>

The shape of the binding curve determines the accuracy of  $K_a$  and  $\Delta H$

# Application: Binding Interactions

**Purpose:** Measure the binding affinity between a ligand and enzyme

## Method:


- (A) 120  $\mu\text{M}$  carboxybenzenesulfonamide (CBS) was titrated into 15  $\mu\text{M}$  Carbonic Anhydrase II
- (B) 30  $\mu\text{M}$  acetazolamide (ACTZ) was titrated into 2.5  $\mu\text{M}$  Carbonic Anhydrase II
- Run in duplicates

## Observations:

- The ITC technology generates a complete thermodynamic profile of the interaction
- Different ligands bind with different strength
- Enthalpy and stoichiometry provide additional information

*ITC is can be used to screen drug candidates by binding constants*

## Protein-Ligand Interactions



The shape of the binding curve determines the accuracy of  $K_a$  and  $\Delta\text{H}$

Cartoon created from 6RFH pdb.org

Cartoon from 3HS4 pdb.org

# Application: Binding Interactions

**Purpose:** Measure the binding affinity between a ligand and protein

## Method:



- Ligand was titrated into protein solution

## Observations:

- The ligand binds two chemically distinct sites on a receptor molecule.
- Binding isotherm for this type of interaction can have two inflections as shown in figure

*The shape of the curve depends on the affinities of the ligand.*

## Protein-Ligand: Multiple Binding Sites



The shape of the binding curve suggests multiple binding sites

# Application: Binding Interactions

**Purpose:** Measure the energetics of dissociation of bovine insulin in aqueous solution

## Method:


- A concentrated solution of protein is injected into buffer in the cell

## Observations:


- Initially, dimer dissociates when it enters the cell, since its concentration in the cell is low.
- As its concentration increases, the dimers dissociate less, until eventually it reaches a concentration that drives dimer formation.
- The dimer dissociation model is used to study dissociation events.  $K_d$  in this model represents the dissociation constant:  $K_d = 22 \mu\text{M}$

*ITC is also suited for studying the dissociation of homodimers*

## Dimer Dissociation



## Human Insulin



Data is consistent with the ITC literature values for bovine insulin collected at pH 2.5

# Application: Binding Interactions

## Purpose: Critical Micelle Concentration (CMC)

### Method:






- Titrate concentrated detergent suspension (micelles) into buffer.

### Observations:

- Critical micelle concentration (CMC) is the concentration at which detergents aggregate to form micelles.
- Initially micelles dissociate in sample cell. At CMC, detergent in the sample cell aggregates. Midpoint of the inflection is the CMC.

*ITC can provide a comprehensive description of surfactant association*

## Critical Micelle Concentration (CMC)



# Application: Binding Interactions

**Purpose:** Evaluate enzyme kinetics under steady-state conditions

## Method:

- Fill syringe with excess substrate, the amount added at the end of the titration will need to be greater than the  $K_m$ .
- Fill cell with a limiting amount of enzyme.

$$[S]_{\text{cell final}} > K_m$$


## Observations:

- Substrate concentration-dependent heat flow of enzymatic reactions may be used for kinetics analysis and determination of Michaelis-Menten reaction parameters
- The plateau is related to the maximum turnover, or  $v_{\text{max}}$ . Additional substrate does not increase turnover at this point, the enzyme is working at its maximum capacity.


*ITC is a facile technique for characterizing enzyme kinetics parameters*

## Enzyme Kinetics: Multiple Injection Method

Rate of heat flow and reaction enthalpy



Michaelis-Menten and Lineweaver-Burk plots



Multiple injections of substrate into enzyme allows for multiple rate determinations under steady-state conditions within one experiment

**Purpose:** Investigate if biotherapeutic is still active?

## Method:

- Titrate known ligand solution into biotherapeutic solution

## Observations:


- Check for unexpected stoichiometry

## Uses:

- Biotherapeutic quality
- This technique is currently used by customers to test antibody batches and the effects of immobilization

*ITC can be a useful tool in evaluating loss of activity in biotherapeutics*

## Anti-quinidine antibody batches compared



Waters™



# NanoDSC Theory and Instrument

# Why use a Nano DSC?

- DSC is the only technique that allows the direct measure of  $T_m$ ,  $DC_p$  and  $DH$ .
  - DSC allows for calculation of entropy ( $DS$ ) and free energy ( $DG$ ).
- Does not require nitrogen purge
  - i.e. CD-spectroscopy
- Equally useful for macromolecules and small molecules (no MW limit)
- Compatible with essentially any buffer or additive
- Requires small sample concentrations and volume



# Nano DSC Specifications and Configurations

|                              |  |
|------------------------------|--|
| <b>Temp Range</b>            | <b>-10° C to 130° C Capillary<br/>-10° C to 160° C Cylindrical</b> |
| <b>Pressure Perturbation</b> | <b>To 6 atm</b>  |
| <b>Cell volume</b>           | <b>0.3 mL</b>  |
| <b>Cell material</b>         | <b>Platinum – Capillary<br/>Gold – Cylindrical</b>                 |

Capillary can also be automated

# Nano DSC Cell Geometry

Continuous Capillary (130°) Fixed-cell for maximum sensitivity



- Cell Construction; Inert to biomaterials 99.99% Platinum
- Small Sample Volume (0.3 mL)
- Attenuates or delays onset of aggregation until after protein has unfolded
- Easy-to-fill and clean design





# Benefit of Capillary Sample Cell

Data obtained with a DSC with “Coin and Cylindrical” shaped sample cell




Data obtained with a Nano DSC with continuous Capillary sample cell



Purified human IgG<sub>1</sub> monoclonal antibody in physiological buffer;  
0.5 mg/mL

# Two state model of protein unfolding



## Enthalpically Favorable

Hydrophobic & electrostatic interactions,  
H<sup>+</sup> bonding

## Entropically Favorable

Changes in solvation &  
conformational freedom

- Heat associated with unfolding (endothermic) and folding (exothermic) is easily measured by calorimetry, allowing thermodynamic analysis of the folding/unfolding process.
- Folding and unfolding of a small protein, a domain, or a subunit, is 'cooperative' (once started, it goes to completion).
- These small units can fold and unfold reversibly. Reversibility is directly measurable using DSC.

## ■ Heat capacity change ( $\Delta C_p$ )


- Trend: Protein with lower  $\Delta C_p$  - more rigid and protein with higher  $\Delta C_p$  - more flexible.
- Primarily reflects exposure of hydrophobic groups.  $\Delta C_p$  is positive.

## ■ Enthalpy

- $\Delta H$ : e<sup>2</sup> non-covalent interactions – hydrophobic & electrostatic, H+ bonding

## ■ $T_m$

- Indicates macromolecular stability
- From these measured values it is possible to calculate the entropy ( $\Delta S$ ) and free energy ( $\Delta G$ )



ASTM E2603-08 - for verification of enthalpy and temperature of a fixed-cell DSC ([www.astm.org](http://www.astm.org))


Waters™



# NanoDSC Applications

- Determine thermal transition (melting) temperatures
- Measure  $\Delta H$  of denaturation
- Measure reversibility of thermal processes
- Measure  $\Delta C_p$  of unfolding
- Determine stability macromolecules
- Measure high affinity binding (up to  $10^{20} \text{ M}^{-1}$ )
- Investigate a complex milieu


- Excipient to stabilize against chemical and physical degradation
- Choice of an additive or a formulation is generally determined empirically
- DSC is the fastest way of evaluating additives effect on  $T_m$ , reversibility



Olsen et al., *Thermochimica Acta*, **484**, (2009), 32-37

Many users stop data analysis after reporting a  $T_{\text{max}}$ , the following slides will discuss what is being missed.


- Mutate proteins: more stable, more specific, faster, new properties, etc.
- Useful to predict outcome of a mutation, so need a database of thermodynamically characterized proteins
- Complicated network of interactions
  - Example: enthalpic changes (changes in hydrophobic interactions, hydrogen bonding, electrostatic interactions) are compensated by entropic changes (changes in solvation, conformational freedom)




Grzesiak et al., J. Mol. Bio. **301**, 205-217, 2000

- Since  $\Delta G = \Delta H - T\Delta S$ , unfolding occurs when  $T\Delta S$  increases sufficiently (e.g. by absorbing heat) to overcome stabilizing enthalpic interactions
- Biopolymer unfolding is endothermic

# How to determine if unfolding is kinetically controlled?




- Scan rate dependence of  $T_m$  indicates that N and U are not in equilibrium.
- Their concentrations change at a rate equal to the sum of the unfolding and refolding reactions. Increasing temperature faster than system responds distorts  $T_m$  and shape.




Remmele et al., *JACS* **127**, (2005), 8328-8339

- Interleukin-1 receptor scanned at 0.25, 0.5, 1.0 and 1.5 °C/min:
  - Scan-rate dependence of  $T_m$  indicates folded and unfolded protein are not in equilibrium
  - Unfolding is kinetically controlled

# How to determine if unfolding is reversible?




- Aggregation is characterized by thermodynamic and kinetic components.
  - Thermodynamic component causes unfolding.
  - Kinetic component can result in partial or complete irreversibility.



- Non-reversibility of protein unfolding indicates:
  - Multi-domain or subunit structure
  - Chemical alterations to the sequence
  - Kinetic events hindering



1<sup>st</sup> heating scan to 6<sup>th</sup> heat scan shown



- Unfolding of domains and subunits with different thermal stabilities may produce asymmetric thermograms.
- A small change in sequence, or other alteration, can affect the stability of the whole protein, or the stability of one domain or subunit.
- DSC quickly reveals these stability changes. Practical implications: Identify subunit stability within a protein.

Medved *et al.*, *J. Biol. Chem.* **274**, 717-727, 1999


- Linearized pBR 322 heated at 0.1 °C/min for optical profile, 1 °C/min for DSC scan.
- Individual spectral peaks reflect the cooperative melting of one or more domains.
- DSC provides a complimentary and more direct means of measuring thermodynamic parameters of DNA.



Volker *et al.*, *Biopolymers* **50**, 303-318, 1999



# Importance of concentration determination

- Effect of sample concentration dependence of  $T_m$  is a test for oligomerization.





# Membrane proteins and membranes: detergents and lipids

- Membrane proteins must be solubilized in detergents or lipids
- Detergents form nanospheres in water composed of a single layer of detergent molecules. This is called a micelle.
- Lipids form nanospheres in water composed of a double layer of lipid molecules, called a liposome.




- Many proteins are associated with membranes and serve numerous functions.
- How does membrane stability change in the presence of various proteins
  - DSC can provide indications
- Vesicles were prepared by sonication bath with 5 mg/mL DPPC and then treated
  - Palmitoylated SNAP-25 (Orange)
  - 1.2  $\mu\text{g}$  SNAP-25 (Red)
  - 2  $\mu\text{g}$  SNAP-25 (Blue)
  - Untreated control (Green)




**Palmitoylated SNAP-25 increase membrane stability, while non-palmitoylated SNAP-25 decreases membrane stability relative to the control**

- If a ligand binds preferentially to a folded protein, the  $T_m$  of the protein will generally increase. The more bound ligand there is, or the tighter it binds, the more  $T_m$  increases.
- Can determine binding constant at  $T_m$ .
- Useful for very tight binding
- DSC is a quick way to determine if two molecules interact.



Binding of 2'-CMP to RNase A  $\pm$  5% DMSO (black)  
 $K_a = 5900 \text{ M}^{-1}$  (-DMSO);  $6900 \text{ M}^{-1}$  (+DMSO) at  $T_m$

- Protein structures are dynamic and constantly fluctuating between partially folding and folded structures
- Ligands can preferentially bind to a partially unfolded conformation
- Ex: Zn<sup>2+</sup> added to Ca<sup>2+</sup>-saturated α-lactalbumin, a two-domain protein
- Enthalpies of the domains in absence of Zn<sup>2+</sup> are consistent with the crystal structure
- Increasing [Zn<sup>2+</sup>] progressively destabilizes the protein, driving equilibrium towards the unfolded state
- DSC is a quick way of determining if a ligand binds preferentially to folded or partially unfolded protein.



# Thank You

---

The World Leader in Thermal Analysis, Rheology, and Microcalorimetry

