














Fluorescent dyes

	Alexa Fluor 350		Alexa Fluor 647		Oregon Green
	Alexa Fluor 405		Alexa Fluor 660		Fluorescein (FITC)
	Alexa Fluor 488		Alexa Fluor 680		Cy3
	Alexa Fluor 532		Alexa Fluor 700		Pacific Orange
	Alexa Fluor 546		Alexa Fluor 750		PE-Cyanine7
	Alexa Fluor 555		Pacific Blue		PerCP-Cyanine5.5
	Alexa Fluor 561		Coumarin		Tetramethylrhodamine (TRITC)
	Alexa Fluor 568		BODIPY FL		Texas Red
	Alexa Fluor 594		Pacific Green		Cy5



eBioscience Super Bright and eFluor Dyes

	eFluor 450		PerCP-eFluor 710		Super Bright 645
	eFluor 506		APC-eFluor 780		Super Bright 702
	eFluor 660		Super Bright 436		Super Bright 780
	PE-eFluor 610		Super Bright 600		

DNA Stains

	DAPI
	SYTOX Green SYTO 9
	TO-PRO-3 Propidium Iodide

Fluorescent protein labels

	R-Phycoerythrin (R-PE)
	Allophycocyanin (APC)

Expressed fluorescent proteins

	CFP
	GFP (emGFP)
	RFP (tagRFP)